

FEATURES:

- Lightweight titanium handles
- Hardened stainless steel front-ends
- Low pressure tip activation
- Single piece spring construction
- Color coded by gauge
- Tip ID on handles
- Consistent grip to tip lengths
- Countersunk traction grip design
- Supplied with a sterilizing case



*Every instrument is supplied
with a sterilizing case and cleaning hub.*

ORDERING INFORMATION:

K5-7703	Forceps, Serrated Jaws	23ga.	(teal)
K5-7713	Forceps, End Gripping	23ga.	(teal)
K5-7715	Forceps, End Gripping	25ga.	(gold)
K5-7723	Forceps, Asymmetrical	23ga.	(teal)
K5-7725	Forceps, Asymmetrical	25ga.	(gold)
K5-7733	Forceps, ILM	23ga.	(teal)
K5-7735	Forceps, ILM	25ga.	(gold)
K4-5423	Scissors, curved	23ga.	(blue)
K4-5425	Scissors, curved	25ga.	(blue)

 **katena**
DESIGNED FOR SIGHT™

 **katena**
DESIGNED FOR SIGHT™

4 Stewart Court • Denville, NJ 07834
(973) 989-1600 • (800) 225-1195 • (973) 989-8175 fax
globe@katena.com • www.katena.com

 **katena**
DESIGNED FOR SIGHT™

Vitreo-Retinal INSTRUMENTS



FUNCTION.
FEEL.
AFFORDABILITY.

It's pretty simple, in VR instrumentation you have to get all three of these right, every time.



Katena is proud to introduce a new line of instruments for vitreo-retinal surgery that will exceed your expectations.

Each reusable instrument is produced with Katena quality and subjected to the most stringent standards for flawless *function*. The fit, *feel* and tactile feedback through the surgeon's hand is superb with attention to every detail—weight, balance, grip, diameter and pressure settings.

Lastly, *affordability*, at a price point between current reusables and disposables, these instruments are truly cost effective.


Contact your local sales representative or distributor for more information.

FORCEPS




Serrated Jaws
3mm long

K5-7703 23ga.




End Gripping
2.75mm long

K5-7713 23ga.
K5-7715 25ga.



Asymmetrical
1.5mm long


K5-7723 23ga.
K5-7725 25ga.



ILM
2.75mm long

K5-7733 23ga.
K5-7735 25ga.

SCISSORS



Scissors, curved
1.25mm long

K4-5423 23ga.
K4-5425 25ga.

 **katena**
DESIGNED FOR SIGHT™