

REFRACTION UNIT

TOMEY REFRACTION UNIT TRU-1000

- Space saving (1.55 m²)
- Easy to install (right or left)
- One-man installation
- Electrical lift chair
- Preconfigured for Tomey equipment
- Timeless and modern design
- Low shipment costs (1 box /euro pallet size, 120 x 80 cm)





QUALITY IN DETAIL

> The Tomey Refraction Unit TRU-1000

If space is an issue for you, the Tomey TRU-1000 refraction unit is the perfect solution.

On only 1.55m² you can use this unit for the most common ophthalmic equipment.

Thanks to the sophisticated manufacturing and design, the TRU-1000 offers best versatile performance while heaving all required features for the daily use.




The rotatable control panel enables the user a comfortable operation of the whole unit and the additional devices.

SPECIFICATIONS

Table	rotating 90°, laterally movable for 2 units
Table height	80 cm
Seat height.	43 cm (lowest position)
Seat elevation	43 - 58 cm
Nominal lifting capacity	max. 120 kg
Projector column height	170 cm
Connectors.	2 x 6V/12V, 2 x RS-232, 1 x slit lamp power connector, 1 x USB 2.0
Reading light	white LED 3.5 V
Power consumption	300 VA
Power supply	230 V
Total weight	approx. 190 kg (without devices)



DIMENSIONS



TOMEY EUROPE
TOMEY GmbH
Am Weichselgarten 19a
91058 Erlangen
Germany

TOMEY ASIA-PACIFIC
TOMEY CORPORATION JAPAN
2-11-33 Noritakeshinmachi
Nishi-ku, Nagoya 451-0051
Japan

Phone (+49) - 9131 - 77710
Fax (+49) - 9131 - 777120
eMail: info@tomey.de

Phone (+81) - 52 - 581 - 5327
Fax (+81) - 52 - 561 - 4735
eMail: intl@tomey.co.jp

